



СТОМАТОЛОГИЧЕСКИЕ СВЕТИЛЬНИКИ

Москва (495) 781 00 36 Хабаровск (4212) 46 00 70
dent@dent.rumail002@dent.ru www.dent.ru


Оплата по курсу ЦБ+3%


Product:	Light source:	Brightness:	Dimensions:	Weight:	Price:	Halo
H1200	LED 100 W	200-3000 lux/1,2 m	1200 x 600 x 45 mm	16 kg	3.360 €	
Options: Flat-screen bracket 320 €		Halo is our very first fixture made with the latest LED technology, providing a unique colour reproduction and low energy consumption.				


▲ Colour temperature: 5500K RA90

Product:	Light source:	Brightness:	Dimensions:	Weight:	Price:	Clair
CL104	LED 104-216 W	3000 lux/1,2 m	1530 x 670 x 59 mm	16,3 kg	2.834 €	
Options: Indirect light 233 € Dynamic light 488 €		Clair is a LED fixture with low energy consumption and minimal heat radiation. Large shadow and glare-free workspace. The optional dynamic light is a world first, on a clinic unit.				


▲ Colour temperature: 3000K RA93 - 6000K RA85


Product:	Light source:	Brightness:	Dimensions:	Weight:	Price:	Brite
BR80	LED 85 W	8000 lux/0,6 m	1200 x 235 x 45 mm	5 kg	785 €	
Options: Dimmer with remote control 153 € Tempered glass 68 €		Brite is designed to suit work benches in medical environments as well as offices. Brite is an efficient light source, shadow- and dazzlefree with no heat radiation.				

Product:	Light source:	Brightness:	Dimensions:	Weight:	Price:	Brite Mid
BR40	LED 43 W	4800 lux/0,6 m	640 x 235 x 45 mm	3 kg	548 €	
Options: Dimmer with remote control 153 € Tempered glass 50 €		Brite Mid offers the same benefits as Brite but is suitable for shorter work surfaces.				


Product:	Light source:	Brightness:	Dimensions:	Weight:	Price:	Brite Duo
BR160	LED 170 W	3850 lux/1,2 m	1200 x 650 x 45 mm	2 x 5 kg	1.876 €	
Options: • INCL. REMOTE CONTROL • INCL. DIMMER		Pairing up two of Brite and you get a brightness spread over an even larger area, that is useful above a dental chair or similar.				

▲ Colour temperature: 3000K RA >90 alt. 6500K RA86


Product:	Tubes:	Brightness:	Dimensions:	Weight:	Price:	Primo U
G66U3ELI	4 x 54 W 2 x 24 W	3900 lux/1,2 m	1416 x 1024 x 55 mm	18 kg	2.345 €	
• INCL. INDIRECT LIGHT • INCL. DIMMER (ON THE MIDDLE PART)		Primo U consists of three long, narrow lamps in a U-shape providing an airier light distribution. It complies with the stringent requirements for cleanliness, being easy to maintain and with easily replaced fluorescent tubes.				
Options: Flat-screen bracket 320 € Full Dimmer 278 €						


Product:	Tubes:	Brightness:	Dimensions:	Weight:	Price:	Denta Hybrid
G6654ELI	4 x 54 W 2 x 24 W	4000 lux/1,2 m	1560 x 780 x 55 mm	22 kg	2.545 €	
Options: LED Light 985 € Dimmer 440 € Flat-screen bracket* 320 €		Denta Hybrid is a pioneering lighting system offering a unique combination of fluorescent tubes and LED lighting that allows colours to change at the speed you select. *Cannot be used in combination with LED Light.				

Colour temperature of tubes: 6500K RA93


Product:	Tubes:	Brightness:	Dimensions:	Weight:	Price:	Denta CDP
T5628EL	6 x 28 W	4000 lux/1,2 m	1355 x 700 x 55 mm	11 kg	1.665 €	
Options: Indirect light (2 x 28 W) 265 €		Denta is our standard model equipped with Conical De-glaring Prism, providing a soft shadow- and dazzle-free light.				

▲ Colour temperature of tubes: 6500K RA85 alt. 4000K RA85

Product:	Tubes:	Brightness:	Dimensions:	Weight:	Price:	Dentaplus CDP
T5654ELD	6 x 54 W	280-5800 lux/1,2 m	1355 x 700 x 55 mm	11 kg	2.160 €	
• INCL. DIMMER Options: Indirect light (2 x 54 W) 295 €		Dentaplus is equipped with a dimmer switch that adjusts the brightness to achieve desired effect.				

Product:	Tubes:	Brightness:	Dimensions:	Weight:	Price:	Denta Mid
T5454EL	4 x 54 W	4500 lux/1,2 m	1355 x 525 x 55 mm	9 kg	1.455 €	
Options: Indirect light 115 € Dimmer 220 €		Denta Mid is a smaller and lighter version of Dentaplus with lighting power reaching up to 4500 Lux.				

▲ Colour temperature of tubes: 6500K RA93 alt. 4000K RA93/RA85

Optional equipment:	Price:	Flat-screen bracket
Halo, Denta Hybrid*, Primo U	320 €	
<p>An easy-to-install flat-screen bracket designed for Halo, Denta Hybrid* and Primo U, which can be added to new or existing lighting systems.</p> <p>*Cannot be used in combination with LED Light.</p>		

Product:	Price:	Mounting equipment/accessories
Wire kit WT5	95 €	Wire cables (1,5 m) incl. 4 ceiling attachments and 4 grippers.
Carrier rail 600	40 €	Two 600 mm carrier rails incl. 4 clips for ceiling mounting.
Carrier rail 1200	55 €	Two 1200 mm carrier rails incl. 6 clips for ceiling mounting.
Wall mounted dimmer POT2	60 €	Push device instead of standard cord dimmer.

Оплата по курсу ЦБ+3%

Prices include light source (LED and/or fluorescent tubes), suspension cables (wires) and packaging (cardboard). Prices are valid ex works Nossebro, excl. VAT and pallets are charged net. All prices quoted are recommended retail prices and reflect the manufacturer's understanding of current market prices. The price list should only be taken as an indication, as each retailer has the right to set his own prices. D-TEC AB reserves the right to make price adjustments. A 5 year warranty is issued on all LED-lighting fixtures (excl. Denta Hybrid LED) and a 2 year warranty on all fixtures with fluorescent tubes.